

CAMARO PP RED

HB-Nr. 499294 • 276000364841869

06.10.2022 aAa: 423 Beta casein: A1A2 Kappa casein: AE ZW: gen.ZW 08.2024 Breeder: Göpel

Pedigree

CARTOON P v. ATOP-RED
SYBILLE P VG-85

LUPO PP v. LUCKY PP
SPARKY GP-84

SPARK RED
SASKIA P GP-84

Total index **RZ€ 1455** **RZG 137**

Production **Si. 73%** **RZM 117**

Milk	+431		
Fat %	+0,13	Protein %	+0,05
Fat kg	+31	Protein kg	+20

Functionality

RZN	Productive Life	121
RZS	Somatic cell score	110
RZR	Daughter Pregnancy Rate	108
RZPersistenz	Persistenz	102
RZD	Milking Speed	107
KVd	Sire calving ease	114
RZK _m	Daughter calving ease	110
RZOeko	Ecological Index	130
RZFE1	Feed Efficiency	108
RZRobot	RZRobot	120

Health

RZGesund 119

RZEuterfit	Udder health	111
RZKlauen	Claw health	107
RZRepro	Fertility	110
RZMetabol	Metabolic stability	106
DDcontrol	DDcontrol	110
RZKälberfit	Calves fitness	114

Exterieur

RZE 117

Dairy Character	96
Body Composite	86
Feet & Legs	112
Udder	124



Camaro PP Red

Camaro PP Red comes from an alternative dam line and is the only red PP sire in Germany at current breeding value level without Solitair and Match blood in the pedigree. In addition, other widely used sires such as McDonald P, Hotspot P, Simon P or Mission P and their sons are not represented.

Camaro convinces with a good overall breeding value with solid performance as well as high longevity and excellent fitness and health values. Camaro inherits less size, good body condition (BCS 107) as well as very good feet & legs and udder. Camaro is the #1 PP sire in Germany after RZRobot (120), can be used on heifers without any problems and is very suitable for organic farms (RZÖko 130).

Linear

Merkmal	ZW	Tendenz	88	112	124	Tendenz
Stature	92	small	<div style="width: 10%;"></div>			tall
Rip structure	101	tight rib	<div style="width: 5%;"></div>			open rib
Body depth	80	shallow	<div style="width: 20%;"></div>			deep
Strength	96	frail	<div style="width: 5%;"></div>			strong
Rump angle	92	high pins	<div style="width: 10%;"></div>			sloped
Thurl width	92	narrow	<div style="width: 10%;"></div>			wide
Rear legs side view	93	posty	<div style="width: 5%;"></div>			sickled
Foot angle	109	flat	<div style="width: 15%;"></div>			steep
Hock quality	96	coarse	<div style="width: 5%;"></div>			dry
Rear legs rear view	113	hocked in	<div style="width: 20%;"></div>			parallel
Front feet orientation	105	hocked in	<div style="width: 10%;"></div>			parallel
Locomotion	109	poor	<div style="width: 15%;"></div>			good
Rear udder height	122	deep	<div style="width: 25%;"></div>			high
Udder cleft	111	weak	<div style="width: 15%;"></div>			strong
Front teat placement	101	wide	<div style="width: 5%;"></div>			close
Rear teat placement	102	wide	<div style="width: 5%;"></div>			close
Fore udder attachment	117	loose	<div style="width: 20%;"></div>			strong
Udder balance	105	ht	<div style="width: 10%;"></div>			vt
Udder depth	118	deep	<div style="width: 20%;"></div>			high
Teat length	98	short	<div style="width: 5%;"></div>			long



GGM: Butlerview RUW Saskia GP-84